

Beatle-02P Series

Digital radiography system



Beatle-02P





Beagle-02P Series

It is the extremely compact and lightweight body, straightforward operation and internal battery power supply of Beagle-02P Series that fits for a variety of examination situations, making it ideal portable solution for mobile veterinary X-ray applications.



Reliable imaging

Powerful and stable output delivers consistently high quality images, suitable for animals' radiography, such as equine X-ray.



Exceptionally lightweight

The highly integrated structure allows the system to be easily transported anywhere and freely positioned to meet the needs of veterinary imaging.



Superior battery performance

Internal battery power brings maximum flexibility, while extended battery life of up to 8 hours guarantees all-day clinical use.



Tripod with durability

- Solid construction combined with three leg-angle stops, the system can be set up on almost any surface.
- Extended vertical travel range with 4-section leg lever lock mechanism enables stable imaging at any chosen height.

Transport with ease

- All necessary components, such as the main X-ray unit, detector, image processing workstation and exposure control, can be all packed into the carrying case.
- The integrated package allows you to move quickly from one location to another.



Workflow with efficiency

- The intelligent image workstation with user-friendly interface allows exposure to be prepared within a few clicks.
- Instant imaging, sophisticated image processing and tailored toolkits helps accurate and rapid diagnosis right next to animals.

